



EIF4G1 Antibody

Product Code	CSB-PA007568ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04637
Immunogen	Recombinant Human Eukaryotic translation initiation factor 4 gamma 1 (EIF4G1), partial (1250-1599AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 4 gamma 1 (eIF-4-gamma 1) (eIF-4G 1) (eIF-4G1) (p220), EIF4G1, EIF4F EIF4G EIF4GI
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Target Names	EIF4G1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.