



flaA Antibody

Product Code	CSB-PA673831ZA01CBV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q46113
Immunogen	Recombinant Campylobacter jejuni flaA protein
Raised In	Rabbit
Species Reactivity	Campylobacter jejuni
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flagellin A, flaA, fla1
Immunogen Species	Campylobacter jejuni
Lead Time	Made-to-order (12-14 weeks)
Target Names	flaA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.