



murG Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA330738ZA01BRJ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P37585 |
| Immunogen | Recombinant Bacillus subtilis murG protein |
| Raised In | Rabbit |
| Species Reactivity | Bacillus subtilis |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | UDP-N-acetylglucosamine--N-acetylmuramyl-(pentapeptide) pyrophosphoryl-undecaprenol N-acetylglucosamine transferase (EC 2.4.1.227) (Undecaprenyl-PP-MurNAc-pentapeptide-UDPGlcNAc GlcNAc transferase), murG |
| Immunogen Species | Bacillus subtilis (strain 168) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | murG |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |