



Toxin CfTX-1 Antibody

Product Code	CSB-PA417651ZA01CWM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7L035
Immunogen	Recombinant Chironex fleckeri Toxin CfTx-1,partial (145-456AA)
Raised In	Rabbit
Species Reactivity	Chironex fleckeri
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toxin CfTX-1 (Toxin 1)
Immunogen Species	Chironex fleckeri (Box jellyfish)
Lead Time	Made-to-order
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.