



Delta-theraphotoxin-Hm1a Antibody

Product Code	CSB-PA351577ZA01HGU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60992
Immunogen	Recombinant Heteroscodra maculata Delta-theraphotoxin-Hm1a (1-35AA)
Raised In	Rabbit
Species Reactivity	Heteroscodra maculata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Delta-theraphotoxin-Hm1a (Delta-TRTX-Hm1a) (Heteroscodratoxin-1) (HmTx1) (Kappa-theraphotoxin-Hm1a) (Kappa-TRTX-Hm1a)
Immunogen Species	Heteroscodra maculata (Togo starburst tarantula?(Togo starburst baboon spider)
Lead Time	Made-to-order
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.