



NS Antibody

Product Code	CSB-PA365963ZA01IJK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03502
Immunogen	Recombinant Influenza B virus NS protein
Raised In	Rabbit
Species Reactivity	Influenza B virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-structural protein 1 (NS1) (NS1A), NS
Immunogen Species	Influenza B virus (strain B/Lee/1940)
Lead Time	Made-to-order (12-14 weeks)
Target Names	NS
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.