




Miraculin Antibody

Product Code	CSB-PA320179ZA01RES
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13087
Immunogen	Recombinant Synsepalum dulcificum Miraculin (30-220AA)
Raised In	Rabbit
Species Reactivity	Synsepalum dulcificum
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Miraculin (MIR)
Immunogen Species	Synsepalum dulcificum (Miracle fruit) (<i>Richadella dulcifica</i>)

Image



Western Blot

Positive WB detected in Recombinant protein
 All lanes: Miraculin antibody at 1:2000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 24.9 kDa
 Observed band size: 30 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.