



M Antibody

Product Code	CSB-PA333055ZA01RJA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25223
Immunogen	Recombinant Rabies virus M protein
Raised In	Rabbit
Species Reactivity	Rabies virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix protein (Phosphoprotein M2), M
Immunogen Species	Rabies virus (strain CVS-11) (RABV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	M
Product Modify	Unmodified