



NP Antibody

Product Code	CSB-PA323576ZA01ZAB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18272
Immunogen	Recombinant Zaire ebolavirus Nucleoprotein (NP), partial (488-739AA)
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (Protein N), NP
Immunogen Species	Zaire ebolavirus (strain Mayinga-76) (ZEBOV) (Zaire Ebola virus)
Lead Time	Made-to-order
Target Names	NP
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.