



HEV1 Antibody

Product Code	CSB-PA365797ZA01HWI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02877
Immunogen	Recombinant Hevea brasiliensis HEV1 protein
Raised In	Rabbit
Species Reactivity	Hevea brasiliensis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pro-hevein (Major hevein) [Cleaved into: Hevein (allergen Hev b 6); Win-like protein], HEV1
Immunogen Species	Hevea brasiliensis (Para rubber tree) (Siphonia brasiliensis)
Lead Time	Made-to-order (12-14 weeks)
Target Names	HEV1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.