




# Outer capsid glycoprotein VP7 Antibody

<b>Product Code</b>	CSB-PA420195ZA01RIU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	A8D8S8
<b>Immunogen</b>	Recombinant Rotavirus A (strain RVA/Cow/Canada/C486/1977/G6P6[1]) (RV-A) Outer capsid glycoprotein VP7 protein (34-309aa)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rotavirus A (strain RVA/Cow/Canada/C486/1977/G6P6[1]) (RV-A)
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Rotavirus A (strain RVA/Cow/Canada/C486/1977/G6P6[1]) (RV-A)
<b>Target Names</b>	Outer capsid glycoprotein VP7 Antibody

## Image



### Western Blot

Positive WB detected in: recombinant protein

All lanes: Outer capsid glycoprotein VP7 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36 kDa

Observed band size: 36 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.