



Lyn Antibody

Product Code	CSB-PA013259ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25911
Immunogen	Recombinant Mus musculus Lyn protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase Lyn (EC 2.7.10.2) (V-src-1 Yamaguchi sarcoma viral related oncogene homolog) (p53Lyn) (p56Lyn), Lyn
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Lyn
Product Modify	Unmodified