



D-2-hydroxyisocaproate dehydrogenase Antibody

Product Code	CSB-PA322950ZA01LAY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17584
Immunogen	Recombinant Lactobacillus paracasei D-2-hydroxyisocaproate dehydrogenase protein
Raised In	Rabbit
Species Reactivity	Lactobacillus paracasei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	D-2-hydroxyisocaproate dehydrogenase (D-HICDH) (EC 1.1.1.-)
Immunogen Species	Lactobacillus paracasei
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.