



ASAH1 Antibody

Product Code	CSB-PA002169ZA01EQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5A6P2
Immunogen	Recombinant Pan troglodytes ASAH1 protein
Raised In	Rabbit
Species Reactivity	Pan troglodytes
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acid ceramidase (AC) (ACDase) (Acid CDase) (EC 3.5.1.23) (Acylsphingosine deacylase) (N-acylsphingosine amidohydrolase), ASAH1
Immunogen Species	Pan troglodytes (Chimpanzee)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ASAH1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.