



# Ccn1 Antibody

<b>Product Code</b>	CSB-PA006463ZA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P18406
<b>Immunogen</b>	Recombinant Mus musculus Ccn1 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mus musculus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein CYR61 (3CH61) (CCN family member 1) (Cysteine-rich angiogenic inducer 61) (Insulin-like growth factor-binding protein 10) (IBP-10) (IGF-binding protein 10) (IGFBP-10), Cyr61, Ccn1 Igfbp10
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	Ccn1
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.