



Pparg Antibody

Product Code	CSB-PA018424ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37238
Immunogen	Recombinant Mouse Peroxisome proliferator-activated receptor gamma (Pparg) (1-505AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisome proliferator-activated receptor gamma (PPAR-gamma) (Nuclear receptor subfamily 1 group C member 3), Pparg, Nr1c3
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order
Gene Names	Pparg
Product Modify	Unmodified