



Probable non-specific lipid-transfer protein 2 Antibody

Product Code	CSB-PA345923ZA01ESB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55958
Immunogen	Recombinant Parietaria judaica Probable non-specific lipid-transfer protein 2 protein
Raised In	Rabbit
Species Reactivity	Parietaria judaica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable non-specific lipid-transfer protein 2 (LTP 2) (Allergen Par j II) (Major pollen allergen Par j 2.0101) (Protein P2) (allergen Par j 2.0101)
Immunogen Species	Parietaria judaica (Pellitory-of-the-wall?(Parietaria diffusa)
Lead Time	Made-to-order (12-14 weeks)
Target Names	LTP2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.