



D1 Antibody

Product Code	CSB-PA814081ZA01NDF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8GTM0
Immunogen	Recombinant Nicotiana alata D1 protein
Raised In	Rabbit
Species Reactivity	Nicotiana alata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flower-specific defensin (NaD1), D1
Immunogen Species	Nicotiana alata (Winged tobacco) (Persian tobacco)
Lead Time	Made-to-order (12-14 weeks)
Target Names	D1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.