



Cyanovirin-N homolog Antibody

Product Code	CSB-PA692568ZA01TJE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5MK11
Immunogen	Recombinant Tuber borchii Cyanovirin-N homolog protein
Raised In	Rabbit
Species Reactivity	Tuber borchii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyanovirin-N homolog (CV-N homolog)
Immunogen Species	Tuber borchii (White truffle)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.