



ystB Antibody

Product Code	CSB-PA303368ZA01YAAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P74977
Immunogen	Recombinant Yersinia enterocolitica ystB protein
Raised In	Rabbit
Species Reactivity	Yersinia enterocolitica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat-stable enterotoxin B (Y-STB), ystB
Immunogen Species	Yersinia enterocolitica
Lead Time	Made-to-order (12-14 weeks)
Gene Names	ystB
Product Modify	Unmodified