



GNAT2 Antibody

Product Code	CSB-PA009599ZA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04696
Immunogen	Recombinant Bos taurus GNAT2 protein
Raised In	Rabbit
Species Reactivity	Bos taurus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein G(t) subunit alpha-2 (Transducin alpha-2 chain), GNAT2
Immunogen Species	Bos taurus (Bovine)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	GNAT2
Product Modify	Unmodified