



LY6G6D Antibody

Product Code	CSB-PA013246ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95868
Immunogen	Recombinant Human Lymphocyte antigen 6 complex locus protein G6d (LY6G6D) (20-104AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte antigen 6 complex locus protein G6d (Protein Ly6-D) (Megakaryocyte-enhanced gene transcript 1 protein), LY6G6D, C6orf23 G6D MEGT1 NG25
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Target Names	LY6G6D
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.