




A27L Antibody

Product Code	CSB-PA324948ZA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20535
Immunogen	Recombinant Vaccinia virus 14 kDa fusion protein (A27L) (1-110AA)
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	14 kDa fusion protein
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Target Names	A27L

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: A27L antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 19.6 kDa

Observed band size: 23 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.