



Thrombin-like enzyme gabonase Antibody

Product Code	CSB-PA314359ZA01BGW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C577
Immunogen	Recombinant Bitis gabonica Thrombin-like enzyme gabonase protein
Raised In	Rabbit
Species Reactivity	Bitis gabonica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombin-like enzyme gabonase (SVTLE) (EC 3.4.21.55) (Fibrinogen-clotting enzyme) (Snake venom serine protease) (SVSP) (Venombin-AB) (Fragment)
Immunogen Species	Bitis gabonica (Gaboon adder) (Gaboon viper)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.