



EBNA1 Antibody

Product Code	CSB-PA365878ZA01EFA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03211
Immunogen	Recombinant Epstein-Barr virus EBNA1 protein
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epstein-Barr nuclear antigen 1 (EBNA-1) (EBV nuclear antigen 1), EBNA1
Immunogen Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Lead Time	Made-to-order (12-14 weeks)
Target Names	EBNA1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.