



ompA Antibody

Product Code	CSB-PA310272ZA01CJB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00087
Immunogen	Recombinant Chlamydomophila psittaci ompA protein
Raised In	Rabbit
Species Reactivity	Chlamydomophila psittaci
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane porin (MOMP), ompA, omp1
Immunogen Species	Chlamydomophila psittaci (Chlamydia psittaci)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ompA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.