



# NCU05495 Antibody

<b>Product Code</b>	CSB-PA742424ZA01NHA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q7S6U4
<b>Immunogen</b>	Recombinant Neurospora crassa NCU05495 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Neurospora crassa
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cyanovirin-N homolog (CV-N homolog)
<b>Immunogen Species</b>	Neurospora crassa (strain ATCC 24698 / 74-OR23-1A / CBS 708.71 / DSM 1257 / FGSC 987)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	NCU05495
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.