



Grem1 Antibody

Product Code	CSB-PA009892ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35793
Immunogen	Recombinant Rat Gremlin-1 (Grem1) (25-184AA)
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gremlin-1 (Cysteine knot superfamily 1, BMP antagonist 1) (Down-regulated in Mos-transformed cells protein), Grem1, Cktsf1b1 Drm
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order
Gene Names	Grem1
Product Modify	Unmodified