



gox Antibody

Product Code	CSB-PA319481ZA01DOZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13006
Immunogen	Recombinant Aspergillus niger gox protein
Raised In	Rabbit
Species Reactivity	Aspergillus niger
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose oxidase (EC 1.1.3.4) (Beta-D-glucose:oxygen 1-oxido-reductase) (Glucose oxyhydrase) (GOD), gox
Immunogen Species	Aspergillus niger
Lead Time	Made-to-order (12-14 weeks)
Target Names	gox
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.