



Ins1 Antibody

Product Code	CSB-PA355622ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01322
Immunogen	Recombinant Rat Insulin-1 (Ins1), partial (25-54AA)
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-1 [Cleaved into: Insulin-1 B chain; Insulin-1 A chain], Ins1, Ins-1
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Ins1
Product Modify	Unmodified