



CLCA1 Antibody

Product Code	CSB-PA886948ZA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9TUB5
Immunogen	Recombinant Sus scrofa CLCA1 protein
Raised In	Rabbit
Species Reactivity	Sus scrofa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-activated chloride channel regulator 1 (EC 3.4.-.-) (Calcium-activated chloride channel family member 1) (pCLCA1), CLCA1, AECC
Immunogen Species	Sus scrofa (Pig)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	CLCA1
Product Modify	Unmodified