



ACTB Antibody

Product Code	CSB-PA344170ZA01DXU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48975
Immunogen	Recombinant Cricetulus griseus ACTB protein
Raised In	Rabbit
Species Reactivity	Cricetulus griseus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, cytoplasmic 1 (Beta-actin) [Cleaved into: Actin, cytoplasmic 1, N-terminally processed], ACTB
Immunogen Species	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ACTB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.