



SNCG Antibody

Product Code	CSB-PA646889ZA01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2PFW6
Immunogen	Recombinant Macaca fascicularis SNCG protein
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-synuclein, SNCG
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Lead Time	Made-to-order (12-14 weeks)
Target Names	SNCG
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.