



Outer capsid glycoprotein VP7 Antibody

Product Code	CSB-PA325524ZA01RFU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17700
Immunogen	Recombinant Rotavirus A Outer capsid glycoprotein VP7 protein
Raised In	Rabbit
Species Reactivity	Rotavirus A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer capsid glycoprotein VP7
Immunogen Species	Rotavirus A (isolate RVA/Cow/United States/B223/1983/G10P8[11]) (RV-A)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.