



Klra3 Antibody

Product Code	CSB-PA720664ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q64329
Immunogen	Recombinant Mus musculus Klra3 protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell lectin-like receptor 3 (5E6) (Lymphocyte antigen 49c) (Ly-49c) (Nk2.1) (T-cell surface glycoprotein Ly-49C), Klra3, Ly-49c Ly49C
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Klra3
Product Modify	Unmodified