



Outer capsid protein VP4 Antibody

Product Code	CSB-PA339583ZA01RFU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35746
Immunogen	Recombinant Rotavirus A Outer capsid protein VP4 protein
Raised In	Rabbit
Species Reactivity	Rotavirus A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer capsid protein VP4 (Hemagglutinin) [Cleaved into: Outer capsid protein VP8*; Outer capsid protein VP5*]
Immunogen Species	Rotavirus A (isolate RVA/Cow/United States/B223/1983/G10P8[11]) (RV-A)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.