




Alpha-conotoxin VxXXC Antibody

Product Code	CSB-PA314942ZA01DWK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C1W7
Immunogen	Recombinant Conus vexillum (Flag cone) Alpha-conotoxin VxXXC protein (1-47aa)
Raised In	Rabbit
Species Reactivity	Conus vexillum (Flag cone)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Conus vexillum (Flag cone)
Target Names	Alpha-conotoxin VxXXC Antibody

Image



Western Blot

Positive WB detected in: recombinant protein
All lanes: Alpha-conotoxin VxXXC Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 22 kDa
Observed band size: 22 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.