



nifH Antibody

Product Code	CSB-PA321828ZA01RKN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20623
Immunogen	Recombinant Rhizobium leguminosarum nifH protein
Raised In	Rabbit
Species Reactivity	Rhizobium leguminosarum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitrogenase iron protein (EC 1.18.6.1) (Nitrogenase Fe protein) (Nitrogenase component II) (Nitrogenase reductase) (Fragment), nifH
Immunogen Species	Rhizobium leguminosarum
Lead Time	Made-to-order (12-14 weeks)
Target Names	nifH
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.