



C1qa Antibody

Product Code	CSB-PA003637ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98086
Immunogen	Recombinant Mouse Complement C1q subcomponent subunit A (C1qa) (23-245AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C1q subcomponent subunit A, C1qa
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order
Gene Names	C1qa
Product Modify	Unmodified