



FOXO3 Antibody

Product Code	CSB-PA008836ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43524
Immunogen	Recombinant Homo sapiens FOXO3 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein O3 (AF6q21 protein) (Forkhead in rhabdomyosarcoma-like 1), FOXO3, FKHL1 FOXO3A
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	FOXO3
Product Modify	Unmodified