



UL128 Antibody

Product Code	CSB-PA321481ZA01HWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16837
Immunogen	Recombinant Human cytomegalovirus UL128 protein
Raised In	Rabbit
Species Reactivity	Human cytomegalovirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein UL128, UL128
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (Human herpesvirus 5)
Lead Time	Made-to-order (12-14 weeks)
Target Names	UL128
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.