



LDBPK_361420 Antibody

Product Code	CSB-PA25309ZA01LOE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E9BTK1
Immunogen	Recombinant Leishmania donovani LDBPK_361420 protein
Raised In	Rabbit
Species Reactivity	Leishmania donovani
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transitional endoplasmic reticulum ATPase, putative (LDBPK_361420)
Immunogen Species	Leishmania donovani (strain BPK282A1)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	LDBPK_361420
Product Modify	Unmodified