



IL15RA Antibody

Product Code	CSB-PA614402ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13261
Immunogen	Recombinant Homo sapiens IL15RA protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-15 receptor subunit alpha (IL-15 receptor subunit alpha) (IL-15R-alpha) (IL-15RA) (CD antigen CD215) [Cleaved into: Soluble interleukin-15 receptor subunit alpha (sIL-15 receptor subunit alpha) (sIL-15R-alpha) (sIL-15RA)], IL15RA
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IL15RA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.