



fno Antibody

Product Code	CSB-PA520438ZA01MSQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D9PVP5
Immunogen	Recombinant Methanothermobacter marburgensis fno protein
Raised In	Rabbit
Species Reactivity	Methanothermobacter marburgensis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	F420-dependent NADP reductase (EC 1.5.1.40) (F420H2:NADP oxidoreductase), fno
Immunogen Species	Methanothermobacter marburgensis (strain ATCC BAA-927 / DSM 2133 / JCM 14651 / NBRC 100331 / OCM 82 / Marburg) (Methanobacterium thermoautotrophicum)
Lead Time	Made-to-order (12-14 weeks)
Target Names	fno
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.