



Non-structural glycoprotein 4 Antibody

Product Code	CSB-PA357389ZA01RIV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08434
Immunogen	Recombinant Rotavirus A Non-structural glycoprotein 4 protein
Raised In	Rabbit
Species Reactivity	Rotavirus A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-structural glycoprotein 4 (NSP4) (NCVP5) (NS28)
Immunogen Species	Rotavirus A (strain RVA/Cow/United States/NCDV-Lincoln/1969/G6P6[1]) (RV-A) (Rotavirus A (strain Nebraska calf diarrhea virus))
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.