



# mitF Antibody

<b>Product Code</b>	CSB-PA304681ZA01NGS
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P67875
<b>Immunogen</b>	Recombinant Neosartorya fumigata mitF protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Neosartorya fumigata
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ribonuclease mitogillin (EC 3.1.27.-) (Allergen Asp f I) (Allergen I/a) (IgE-binding ribotoxin) (Major allergen Asp f 1) (allergen Asp f 1), mitF, aspF1
<b>Immunogen Species</b>	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	mitF
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.