



Zinc finger CCCH domain-containing protein ZFN-like Antibody

Product Code	CSB-PA866006ZA01EWE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SWF9
Immunogen	Recombinant Pisum sativum Zinc finger CCCH domain-containing protein ZFN-like protein
Raised In	Rabbit
Species Reactivity	Pisum sativum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger CCCH domain-containing protein ZFN-like
Immunogen Species	Pisum sativum (Garden pea)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.