



# Milr1 Antibody

<b>Product Code</b>	CSB-PA662302ZA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q3TB92
<b>Immunogen</b>	Recombinant Mus musculus Milr1 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mus musculus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Allergin-1 (Allergy inhibitory receptor 1) (Mast cell antigen 32) (MCA-32) (Mast cell Ag-32) (Mast cell immunoglobulin-like receptor 1), Milr1, Gm885 Mca32
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	Milr1
<b>Product Modify</b>	Unmodified