



alr2278 Antibody

Product Code	CSB-PA24909ZA01FUZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8YUQ7
Immunogen	Recombinant Nostoc sp. alr2278 protein
Raised In	Rabbit
Species Reactivity	Nostoc sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alr2278 protein (Alr2278)
Immunogen Species	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Lead Time	Made-to-order (12-14 weeks)
Target Names	alr2278
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.