



alr2278 Antibody

Product Code	CSB-PA24909ZA01FUZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8YUQ7
Immunogen	Recombinant Nostoc sp. Alr2278 protein (alr2278) (1-189AA)
Raised In	Rabbit
Species Reactivity	Nostoc sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Alr2278 protein (Alr2278)
Species	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Lead Time	Made-to-order
Gene Names	alr2278
Product Modify	Unmodified